COLUMN USER GUIDE for Agilent Reversed-Phase Columns

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Poroshell 120/300 / Pursuit / Polaris / HC-C18(2) / TC-C18(2)

Agilent Technologies

This booklet provides general information for all ZORBAX, Poroshell, Pursuit, and Polaris reversed-phase columns.

For additional detailed information about your specific phase or family, see: **[agilent.com/chem/columnchoices](https://www.agilent.com/chem/columnchoices)**

Getting Started

A QC Column Performance Report, including a test chromatogram, is enclosed with every Agilent column. The QC test system has been modified from a standard system to minimize system dead volume, so it may vary from the system used in your lab. This allows a better evaluation of the column and assures a more consistent product. A properly configured LC system will generate similar results to the chromatogram on your QC Performance Report.

Modern columns are robust and are designed to operate for long periods under normal chromatographic conditions. You can maximize column performance by running it within specifications. Always review the specifications before putting in place a final method.

Using Your Column

Installation

- The direction of flow is marked on the column.
- 1.8 μm columns (ZORBAX RRHT, ZORBAX RRHD) can only be operated in the direction flow marked on the column.
- Agilent recommends using Polyketone fittings (5042-8957) for columns up to 600 bar, and A-Line Quick Connect assembly (5067-5961 for 0.075 x 105 mm) or Quick Turn Fitting (5067-5966) for columns that will be operated at UHPLC pressures. Be sure to order the right capillary connection for your fitting.

Learn more at www.agilent.com/chem/a-linefittings.

*Polyketone fitting,
p/n 5042-8957*

*A-Line Quick Connect
assembly, p/n 5067-5961*

*Quick Turn fitting,
p/n 5067-5966*

Column conditioning

Every column is tested before shipment and shipped in the test eluent. So, for first use, rinsing with water is not necessary. If mobile phase additives are used (such as buffers or ion-pair reagents) it is advisable to do an intermediate flush with a mobile phase of the correct composition, but without these additions. Flushing with 10 to 20 column volumes should help in transitioning to your mobile phase. For shorter chain chemistries (C8, Phenyl, CN, for example), care should be taken to make sure the column has been properly equilibrated prior to use. This will ensure reproducibility and help prevent retention time drifting.

Important safety considerations

- All points of connection in liquid chromatographic systems are potential sources of leaks. Users should be aware of the toxicity or flammability of their mobile phases.
- Because of the small particle size, dry column packings are respirable. Columns should only be opened in a well-ventilated area.
- Please adhere to operating pressure limits noted for each column (see chart). Exceeding these limits will compromise chromatographic performance and could be unsafe.

Other operating tips

- While generally not harmful to the column, reverse flow should be avoided except to attempt removal of clogged frit (see "column care").
- Always use high purity reagents and chromatography grade solvent to prepare your mobile phase. Degas and filter all mobile phase prior to use.
- Disassembling a column will degrade column performance.
- New columns contain a mixture of organic solvents and water. See your QC Performance Report for the solvent composition in your column. Initially, care should be taken not to pass any mobile phase through the column that may cause a precipitate to form.
- Agilent reversed-phase columns are compatible with water and all common organic solvents.
- The use of a guard column is recommended to protect your column and increase its lifetime.
- Columns should not be maintained at elevated pH or elevated temperature when not in use.
- Avoid use of this column outside of recommended pH ranges for column phase (see next page). Expect reduced lifetime when operating outside the recommended pH and temperature ranges.

*If you are using Poroshell 120 or ZORBAX RRHT or RRHD columns, 1.8 μ m, consider using a **Fast Guard for UHPLC** to protect your analytical column.*

See [agilent.com/chem/fastguards](https://www.agilent.com/chem/fastguards) for more information.

Column Operating Parameters: pH and Temperature

Phase	Recommended pH range	Maximum operating temperature
ZORBAX Eclipse Plus C18 and C8, Eclipse XDB-C8, XDB-C18, XDB-Phenyl, Poroshell 120 EC-C18 and EC-C8, ZORBAX Bonus-RP, HC-C18(2) and TC-C18(2)	pH 2.0 to 9.0	60 °C
Poroshell HPH	pH 2.0 to 11.0	60 °C
ZORBAX Extend C18	pH 2.0 to 11.5	60 °C (40 °C at high pH)
ZORBAX ODS, ZORBAX Phenyl, ZORBAX Eclipse Plus Phenyl-Hexyl, Poroshell 120 PFP, Poroshell 120 Phenyl-Hexyl, ZORBAX XDB-CN, ZORBAX Eclipse PAH, All Pursuit and Polaris bonded phases.	pH 2.0 to 8.0	60 °C
ZORBAX Eclipse PAH	pH 2.0 to 8.0	60 °C at pH <6.0
Poroshell 120 SB-C18, StableBond C18, ZORBAX Rx-C8, ZORBAX SB-Aq, ZORBAX SB-Phenyl, SB-C3, Poroshell 300, 300SB-C3	pH 1.0 to 8.0	90 °C (StableBond C18 and Poroshell 120 SB-C18) 80 °C (Rx-C8 and SB-Aq, SB-C3 and SB-Phenyl) 70 °C for pH <5.0; 40 °C for pH 5.0 to 8.0 (Poroshell 300 and 300SB-C3)
ZORBAX TMS	pH 2.0 to 7.0	60 °C

Note: All silica-based packings have some solubility in pH >6 aqueous mobile phases. When using silica-based columns at pH >6, best column lifetime is obtained at lower temperatures (40 °C max) using low buffer concentrations in the range of 0.01 to 0.02 M. Operating at extreme ends of pH and temperature ranges will have a significant impact on column lifetime.

Maximum Operating Pressures – Columns up to 9.4 mm id

Column type	Particle size	Pressure limit
ZORBAX Rapid Resolution	3.5 μm	400 bar (6000 psi)
ZORBAX Rapid Resolution High Throughput (RRHT)	1.8 μm	600 bar (9000 psi)
ZORBAX Rapid Resolution High Definition (RRHD)	1.8 μm	1200 bar (17000 psi)
Poroshell 120, Poroshell HPH	2.7 μm , 4 μm (Poroshell 120)	600 bar (9000 psi)
Poroshell 300	5 μm	400 bar (6000 psi)
All ZORBAX, Pursuit, and Polaris columns, not noted as RRHT or RRHD, and HC-C18(2), TC-C18(2)	3.0 μm , 3.5 μm , 5 μm	400 bar (6000 psi)

Mobile Phase Selection and Operating Temperatures

The bonded stationary phase is nonpolar in nature and is best used with polar mobile phases such as methanol/water or acetonitrile/water mixtures. Increasing the amount of organic component typically reduces the retention time of the sample.

Recommended Starting Gradients

Stationary phase	Usage notes
Most reversed-phase columns	5% methanol or acetonitrile initially and 100% methanol or acetonitrile as the final solvent.
ZORBAX Eclipse PAH*	30 or 40% acetonitrile initially, to 100% acetonitrile as final solvent. It may be necessary to cool column to 15 to 20 °C for improved resolution.
ZORBAX Bonus-RP	Lower concentration of organic mobile phase modifier is required for compound elution, compared to traditional long-chain alkyl stationary phases.
Large molecule columns (e.g., Poroshell 300, ZORBAX 300SB columns)	Separations of polypeptides, proteins, DNA fragments, etc. often use trifluoroacetic or formic acid modifier for pH control and/or as an ion-pairing additive for desired retention and selectivity. Organic buffers such as TRIS are especially useful for intermediate pH applications.

*Best column lifetime is achieved at < 40 °C.

Column Care

The inlet frit on columns with a particle size of 2.7 μm or greater has a nominal porosity of 2 μm . Samples that contain particulate matter will plug the column inlet frit. ZORBAX and Fast Guard columns and hardware kits (as required) are recommended for use with such samples, and are generally recommended for all column use.

See [agilent.com/chem/fastguards](https://www.agilent.com/chem/fastguards) for more information about guard columns.

Cleaning Your Column/Extending Column Life

For columns that can be backflushed (Poroshell 120, ZORBAX columns with particles $>1.8 \mu\text{m}$ and all Pursuit and Polaris columns), start with a stronger (less polar) solvent.

1. Disconnect column from detector and run wash solvents into a beaker.
2. Start with your mobile phase without buffer salts (water/organic).
Run at 10 to 20 column volumes through.
3. Next, use 100% organic (methanol or acetonitrile).
4. Check pressure to see if it has returned to normal; if not, then,
5. Discard column or consider stronger conditions,
e.g., 75% acetonitrile/25% isopropanol
6. Increase to 100% isopropanol, 100% methylene chloride or 100% hexane (if you use methylene chloride or hexane, you will need to flush the column with isopropanol prior to use and before returning to your reversed-phase mobile phase).

For columns with particles $< 1.8 \mu\text{m}$, do not backflush the column - replace the column.

Storage Recommendations

Long-term storage of silica-based, bonded phase columns should be in a pure organic solvent. If the column has previously been used with a buffered mobile phase, the buffer should first be removed by purging the column with 20 to 30 column volumes of a 50:50 mixture of methanol or acetonitrile and water, followed by 20 to 30 column volumes of the pure solvent. Before storing, end-fittings should be tightly capped with end-plugs to prevent packing from drying out. Columns may be safely stored for short periods in most mobile phases. To protect equipment, it is desirable to remove salts from the instrument and column by purging the column with the same mobile phase without the buffer (e.g. using 60:40 ACN/H₂O to remove a 60:40 ACN/0.02 M phosphate buffered mobile phase). Re-equilibration is rapid with the original mobile phase when using this approach and any danger of corrosion from the salts is eliminated.

Tips for Getting the Best Chromatographic Results

- Optimize your instrument by minimizing tubing lengths between components, to reduce extra column volume and band broadening. Use 0.12 mm id red tubing for Fast LC/high efficiency columns. Learn about capillary options at [agilent.com/chem/lccapillaries](https://www.agilent.com/chem/lccapillaries)
- Ensure the data collection rate is optimized for your column. Use a higher collection rate for Fast LC columns (Poroshell 120, RRHT, RRHD).
- Use sample filtration or other sample prep as appropriate for your sample. Learn more at [agilent.com/chem/sampleprep](https://www.agilent.com/chem/sampleprep)
- Use Agilent certified lamps in your LC instruments for best performance.

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Further Information

For help choosing the right column, request the LC Column Selection Guide or the BioHPLC Column Selection Guide

agilent.com/chem/getguides

The LC Column and Sample Prep Navigator is a web-based tool that can help you find the right column for your method, based on your input. Consult **agilent.com/chem/navigator**

Agilent Ordering Information

For more information on our products and services,
visit our web site at **agilent.com**

For technical support and local information,
visit **agilent.com/chem/columnsupport**

To place an order, visit **agilent.com/chem/wheretobuy**

Agilent offers a complete line of sample preparation products
to support LC and LC/MS applications.

The Agilent Bond Elut SPE and Captiva Filtration Sample Prep
family of products offer the widest range of solutions for
every level of sample cleanliness to help you increase
throughput and enhance the quality of your data.

Learn more at **agilent.com/chem/sampleprep**

This information is subject to change without notice.

© Agilent Technologies, Inc. 2014

Printed in Canada, October 6, 2014

820000-999